DK: Brugsanvisning for repositionerbar probefikseringsenhed.

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[Fremstillingsdato] Se emballagen.

WARNINGS:

1. In accordance with applicable laws and regulations, the product is intended for use by trained clinical personnel only.
2. The product is intended for single use only! Do NOT reuse, reprocess or resterilize the device. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or cause functional failure, which may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of cross-infection, including but not limited to the transmission of infectious diseases from one patient to another. Contamination of the device may result in patient injury, illness or death.
3. Do not use this device for any purpose other than the stated indications.
4. As the device contains polycarbonate/acrylonitrile butadiene styrene (PC/ABS), thermoplastic elastomer (TPE), polyamide 6/6.6 and polyamide 11, it must not be used on patients with known hypersensitivity to one or more of these materials, as allergic reactions may occur.

CAUTION:

1. These instructions for use contain essential information for the safe and effective use of this device. Read these instructions for use thoroughly before use, together with the instructions for use for all equipment to be used during the procedure, and use the device as described. Keep these and all associated instructions for use in a safe and readily accessible place. If you have any questions or comments regarding the information in these instructions for use, please contact EchoVice.
2. The patient shall be informed of and give consent to the operation/procedure and the potential risks and complications, which may result in patient injury, illness or death.
3. Read the instructions for use before use.
4. Do not use this device with probes that have dimensions other than the stated compatible dimensions.
5. Do not use this device with bite blocks that have dimensions other than the stated compatible dimensions.

WARRANTY:

EchoVice warrants that reasonable care and prudence have been exercised in the design and manufacture of this device. This warranty replaces and excludes all other warranties not expressly set out herein, whether express or implied by law or otherwise, including but not limited to any implied warranty of merchantability or fitness for a particular purpose.

Handling, storage and cleaning of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond EchoVice's control, directly affect the device and the results obtained from its use. EchoVice's obligation under this warranty is limited to replacement of this device, and EchoVice is not liable for any indirect or consequential loss, damage or expense arising directly or indirectly from the use of this device.

EchoVice does not assume, and does not authorize any other person to assume on EchoVice's behalf, any other or additional liability in connection with equipment that has been reused, reprocessed or resterilized. EchoVice gives no warranties, whether express or implied, including but not limited to warranties of merchantability or fitness for the intended use, for such equipment.

NAME

EchoM

INTENDED USE

The intended use of EchoM is to stabilize transesophageal (TEE) ultrasound probes used during examinations, treatments and operations. EchoM is a hardware solution that is integrated into the patient's bite block.

The device prevents movement of transesophageal ultrasound probes along the longitudinal axis but retains the rotational movement. The device uses a control unit to lock clamps around the probe, thereby fixating the device. The device is intended for use by trained clinical personnel.

STRUCTURE

EchoM consists mainly of an end cap, fixation unit, tube, clips and trigger unit (see figure 1).


Figure 1.

1. End cap, 2. Fixation unit, 3. Tube, 4. Clips, 5. Trigger unit

SPECIFICATION

Compatible diameter for TEE ultrasound probe: Ø11.5 mm ± 1.5 mm

INDICATION

Examinations, treatments and operations involving TEE ultrasound probes.

CONTRAINDICATION

1. Use without a bite block. The device is intended to be inserted into a bite block and must not be placed directly in the patient's mouth.
2. Improper dental alignment, dental implants or other objects in the patient's mouth may make correct placement of the bite block difficult.

POTENTIAL COMPLICATIONS

1. Damage to the patient's teeth or oral cavity.
2. Complications that are not currently known or observed may occur.

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device shall be reported to EchoVice and to the competent authority of the Member State in which the user and/or patient is established.

CLINICAL BENEFIT

EchoM has no direct diagnostic or therapeutic effect. The clinical benefit of EchoM derives from the device's ability to stabilize the TEE ultrasound probe during examinations and procedures. By reducing unintended movements of the probe while maintaining controlled rotational function, the device supports consistent probe placement and reproducible imaging conditions. This helps facilitate the acquisition of clinically relevant images required for diagnostic assessment and decision-making during procedures, thereby supporting the safe and effective use of transesophageal echocardiography.

INSPECTION BEFORE USE AND PREPARATION

1. Before use, check whether the packaging is damaged.
2. After checking the expiry date, unpack the device and remove it from the packaging.
3. Do not use the device if functional deviations or other concerns are identified.
4. The product must not be used in the event of a fault or damage to the packaging or a functional fault. If in doubt, contact EchoVice.

INSTRUCTIONS FOR USE

Compatibility with equipment

The device is compatible with TEE ultrasound probes with an insertion-tube diameter of $\varnothing 11.5 \text{ mm} \pm 1.5 \text{ mm}$. Before use, check that the probe used meets this specification.

The compatibility with the bite block used shall be verified before use. EchoM must be securely placed in the bite block and must not be able to move or become displaced. If there is any doubt about the compatibility between EchoM and the selected bite block, the device must not be used. Instead, contact EchoVice for guidance before continuing use.

Preparation/unpacking

1. Take EchoM out of the packaging.
2. Carry out a visual inspection of all units for damage or defects before use.


Figure 2.

Assembly

1. Attach the trigger unit to the handle of the TEE ultrasound probe, keeping a distance from other parts of the TEE probe (see figure 3).


Figure 3.

2. Attach the tube to the TEE ultrasound probe's insertion tube using the clips (see figure 4).


Figure 4.

3. Remove the end cap from the fixation unit and insert the TEE ultrasound probe into the bite block. Then place the fixation unit correctly in the bite block (see figure 5).


Figure 5.

4. Insert the TEE ultrasound probe to the desired depth (see figure 6).


Figure 6.

5. Reposition the end cap in the fixation unit as shown in the figure (see figure 7).


Figure 7.

Fixation and use

1. Pull the trigger until the TEE ultrasound probe is fixated. A distinct click is heard, indicating that fixation has been achieved (see figure 8). Ensure that the probe is securely fixated before continuing the procedure.


Figure 8.

2. To release the TEE ultrasound probe, pull the trigger until a click is heard (see figure 9).


Figure 9.

3. Manipulate the probe until the desired position and the desired image are achieved.
4. Pull the trigger again to reactivate fixation in the desired position. A click confirms that fixation is activated.
5. Repeat release and fixation as needed during the procedure.

Removal/Disposal

1. Ensure that the TEE ultrasound probe is not fixated before the device is removed.
2. Remove the device from the TEE ultrasound probe and the bite block.

NOTE

1. The product must be used by trained clinical personnel only.
2. Check the expiry date before use. Do not use the device if the expiry date has passed.
3. The product is intended for single use only. Do not reuse the device.
4. Do not use the device on pediatric patients.
5. Do not operate the fixation unit, the tube or the trigger unit with excessive force, as this may cause damage to the device.
6. Do not test the fixation while the TEE ultrasound probe is being inserted. Testing of the fixation must be carried out before or after the probe has been inserted.
7. Check that the TEE ultrasound probe is correctly fixated before use.

STORAGE

Store the device dry, protected from sunlight and at ambient temperature. Do not expose the packaging to organic solvents, ionizing radiation or ultraviolet radiation. Check the expiry date on the packaging label before use. Store in a hygienically correct manner.

DISPOSAL

After use, the device must be handled and disposed of in accordance with applicable medical practice and relevant local, regional and national laws and regulations. Dispose of the device and packaging in accordance with the hospital's guidelines.

SYMBOLS

	Do not reuse		Non-sterile
	Use by		Manufacturer
	Lot number		Date of manufacture
	Caution/dispose of after use		CE-marking
	Serial number		Medical device
	Unique Device Identifier (UDI)		Damaged packaging
	Keep dry		Consult the instructions for use

[PACKAGING] Packed in PETG trays with TYVEK seal.

[Date of manufacture] See the packaging.

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